

MASS SPECTROMETRY REPORT



Sample Description

Analyzed date:

Analyst: SHAO

Sample: M31233- ω -Agatoxin IVA-03

M.W.: 5202.3

Lot. No.:

Instrument

Probe:

Nebulizer Gas Flow:

CDL:

CDL Temp.:

Block Temp.:

Waters ZQ2000

ESI

10 L/min

52.0v

350 °C

120 °C

Probe Bias:

+3.5kv

Detector:

0.05kv

T. Flow:

0.5ml/min

B. Conc.:

50%H₂O/50%ACN