

M6006-Tetracaine-01-LCMS

Sample Information

LCMS


Acquired by : System Administrator
Sample Name : M6006-Tetracaine-01
Method File : B 5%-95%-P-new - 2.6min.lcm

Date Acquired : 2021/7/27 13:15:54
Injection Volume : 10
Vial : 31


Column : Waters SunFire C18 50*4.6mm 5um 2.000ml/min 2.6min
Pump A : 0.1% FA and 10% ACN in H2O
Pump B : 0.1% FA and 10% H2O in ACN

Chromatogram


mAU



mAU



段号1



Integration Results for DAD

PDA Ch1 254nm


Peak#	Ret. Time	Height	Area	Area%
1	1.213	418775	1789508	99.097
2	2.013	9763	16310	0.903
总计		428538	1805818	100.000

PDA Ch2 214nm


Peak#	Ret. Time	Height	Area	Area%
1	1.213	1957810	8802089	98.512
2	2.303	64499	132979	1.488
总计		2022309	8935068	100.000

Mass Spectrum


Rt=1.187min 负




Rt=1.204min 正




Rt=1.564min 正




Rt=1.744min 正




Rt=1.772min 负




Rt=1.879min 正




Rt=1.930min 负



Rt=2.194min 正



Rt=2.245min 负



Rt=2.402min 负

